
# DIGITAL MICROSCOPE

318-C

## PRODUCT DESCRIPTION

Thank you for purchasing this equipment, it is complex to use, so please read this manual in detail to understand how to operate the microscope before using it.


### IMPORTANT NOTES

1. Please fully charge the product before first use. Please do not charge directly through the PC, but choose a 5V-2A adapter.
2. In PC mode, this device only supports USB3.0 or above interface, not USB2.0 interface.
3. This product is not a biological microscope and is not suitable for cell biology research.
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.
5. This device cannot accurately read microscope magnification. It is a combination of digital and optical magnification microscope, the specific magnification effect is subject to the actual picture taken.
6. When the battery is low, both the storage mode and PC camera mode cannot be used normally. This is not a product failure.
7. Do not touch the lens or other optical parts with your hands, otherwise it will cause image blur and affect imaging quality.
8. Do not disassemble the product or component structure by yourself to avoid abnormal equipment use.

### BATTERY INSTRUCTIONS

1. Please use a safe and compliant 5V/2A home charger to charge this device.
2. This device supports charging while in use.
3. If the device is left unused for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-discharging the battery.

## FUNCTION INTRODUCTION

### RIGHT SELECTION KEY

Scroll down

## MODE SWITCH KEY

Short press to switch photo/video, file browsing mode / long press to switch different color display modes

### LED LIGHT SWITCH

Left for off, right for on

### RESET BUTTON

When the product crashes, press to force shutdown

### TF CARD SLOT

Support expansion 0-64G TF card

### CHARGING INDICATOR

When charging the red light is on, when it's done the green light is on

### CHARGING INTERFACE

Input 5V=2A

### LIGHT POWER SUPPLY PORT

Output 5V=1A

### LIGHT CHARGING PORT

Input 5V=1A

### BRIGHTNESS ADJUSTMENT

Left for dim, right for bright

### RIGHT SELECTION KEY

Confirm selection / screen flip


### MENU KEY

Set the relevant parameters

### OK KEY

Confirm selection / screen flip

## WORKBENCH INTRODUCTION


## PRODUCT SPECIFICATIONS

Lens pixels	12.0MP
Magnification	1600X
Focusing mode	Manual
Fixed focus distance	90mm
Video resolution	1920*1080P
Photo resolution	4032*3024P
Video format	MOV
Photo format	JPG

## OTHER

Screen size	10.1" IPS
Screen resolution	1024*600P
Operating systems	Windows XP/ Vista/7/8/10/11, MacOS X 11.0 or higher
TF capacity	Support expansion 0-64GB
Input voltage & current	5V/2A
Battery capacity	3000mAh Lithium battery
Operating time	1 hrs
Charging time	3 hrs
Operating temperature	-10°C ~ 50°C
Camera working temperature	-15°C ~ 75°C
Light source	30 high brightness LEDs (brightness adjustable)
Languages	Simplified Chinese / Traditional Chinese / Japanese / English / Korean / French / German / Spanish / Italian / Russian / Polish

## STEPS TO USE


### INSERT THE TF CARD

Please insert the TF card when power off.

Please be careful not to mix the insertion direction and the front and back of the TF card.

Do not use excessive force when inserting or pulling out the TF card, otherwise the structure of the card slot may be damaged.

## FUNCTION BUTTON

### SETTING PAGE OPERATION INSTRUCTIONS

1. Double click the "MENU" key to enter the setting page.


2. Click "←" "→" to select the required menu item.

3. Select the required item, click the "OK" key to enter, select "←" "→" the desired parameters with and then click the "OK" key to confirm.

4. When finish setting, Click the "←" "→" key to exit the setting mode.

5. The first time using the device, you need to go to "MENU" to select the required language (the factory setting of the device is: English).

## FILE BROWSING MODE OPERATION INSTRUCTIONS

### FILE BROWSING MODE OPERATION INSTRUCTIONS

1. Click "←" twice to enter the file browsing mode, you can view all the saved photos and videos. If you need to delete the current photo or video file, click "MENU" key to enter the delete page.

2. If the current file is in video format, click "OK" key for play/pause function, click left and right key for fast forward and fast rewind function.

## MODE INTRODUCTION

The product has three mode functions, which are as follows:  
 1. Memory mode: connect to PC, read, write and format TF card.  
 2. PC camera mode: Connect to a computer and use it through computer.  
 3. Normal mode: normal use mode of the device.

### MEMORY MODE

Connect the microscope to the computer via USB data cable, select memory mode, then click the OK button to view saved photos and videos on the computer.

### PC CAMERA MODE

Only supports USB3.0 or above interface, does not support USB2.0 interface

Note: When using PC mode.

1. some function keys will not be available.
2. In PC camera mode, it is normal that the device monitor will not display images.
3. When connected, if the software has no image, please repeat plugging and unplugging the Type-c data cable or replace the cable.

## WINDOWS SYSTEM

Support Windows vista/XP/7/8/10/11 or later.  
**I. Software Installation**  
 1. For WIN 7/8 systems, download and install the "AMCAP" software.  
 2. For IN 10/11 systems, search directly for the software "Windows Camera" (Note: Please disable the default laptop camera in Windows! And you must change the camera's privacy settings, which is needed to allow access).  
 3. Or download the application "Smart Camera" for your computer, download link: [www.inskam.com/download/camera.zip](http://www.inskam.com/download/camera.zip)  
 4. Connect your microscope to the USB port on the back of your computer case with the Type-C cable that comes with the microscope, press and hold the power button to turn on your computer, select to enter PC camera mode and open the computer software to use.  
 5. Different Windows systems are recommended to use different software, please choose to download the corresponding software to use according to the system.


### MACOS SYSTEM

Support MacOS: MacOS X 11.0 or higher.  
**II. the software page display:**  
 1. Software "AMCAP" page, as follows.  
 2. "Windows Camera" page, as follows.  
 3. Software "Smart Camera" page, as follows.

## FAQ HELP

1. **The product cannot be switched on.**  
 Confirm whether the power is connected.  
 Confirm that the power supply is not sufficient due to the computer being used for charging.
2. **The product is dead.**  
 Press the device reset button to force shutdown.
3. **The image is not clear.**  
 Whether the object being measured is correctly positioned directly below the lens.  
 Whether the focal length is adjusted to the optimum position?  
 Whether the light is used correctly?
4. **Product charging.**  
 Please use 5V-2A adapter to charge the device.
5. **TF card cannot be recognized.**  
 Please make sure that the TF card format is FAT32. If not, please format it. If the device does not recognize the card, please reinsert it. The device only supports high-speed storage cards of Class 10 or above, not including Class 10.
6. **When using PC camera mode, the software does not show the picture after connecting to the computer.**  
 Reconnect the product to the computer.  
 Check if you clicked into PC mode.  
 Whether the battery is low.  
 Make sure that you are using the correct data cable and that you are not using a charging cable without data transfer function.

## COPYING PICTURE AND VIDEO FILES

### CHARGING GUIDE

Connect the device to the Type-C / USB adapter for charging (5V=2A), fast charging is not supported.

During charging, the indicator shows a red light.

When the device is fully charged, the indicator shows a green light.

When the battery is low, both memory mode and PC camera mode cannot be used normally, and the battery must be charged beyond the 1/3 position before it can be used again.

## WEEE COMPLIANCE

Proper disposal of this product (Waste Electrical and Electronic Equipment) (for countries with independent collection systems)  
 This marking on the product, accessories or literature indicates that the product and its electronic accessories should not be disposed of with other household waste at the end of their useful life.  
 To prevent uncontrolled waste disposal that may be hazardous to the environment or human health, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources.  
 Household users should contact the retailer or local government office where they purchased this product for details of where they can take these items and how to recycle them in an environmentally friendly manner.  
 Business users should contact their suppliers and check the terms and conditions of their purchase contracts.  
 This product and its electronic accessories should not be mixed with other commercial waste for disposal.

## FCC STATEMENT

Any changes or modifications not expressly approved by the party responsible for compliance may void the user's authority to operate the equipment.  
 This equipment complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
 (1) This equipment will not cause harmful interference.  
 (2) This equipment must accept any interference received, including interference that may cause unintended operation.  
 FCC Radiation Exposure Statement  
 This equipment complies with the FCC radiation exposure limits established for uncontrolled environments.

### WEEE COMPLIANCE

Proper disposal of this product (Waste Electrical and Electronic Equipment) (for countries with independent collection systems)

This marking on the product, accessories or literature indicates that the product and its electronic accessories should not be disposed of with other household waste at the end of their useful life.

To prevent uncontrolled waste disposal that may be hazardous to the environment or human health, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources.

Household users should contact the retailer or local government office where they purchased this product for details of where they can take these items and how to recycle them in an environmentally friendly manner.

Business users should contact their suppliers and check the terms and conditions of their purchase contracts.

This product and its electronic accessories should not be mixed with other commercial waste for disposal.

## PRODUCT DESCRIPTION

Thank you for purchasing this equipment, it is complex to use, so please read this manual in detail to understand how to operate the microscope before using it.


### IMPORTANT NOTES

1. Please fully charge the product before first use. Please do not charge directly through the PC, but choose a 5V-2A adapter.
2. In PC mode, this device only supports USB3.0 or above interface, not USB2.0 interface.
3. This product is not a biological microscope and is not suitable for cell biology research.
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.
5. This device cannot accurately read microscope magnification. It is a combination of digital and optical magnification microscope, the specific magnification effect is subject to the actual picture taken.
6. When the battery is low, both the storage mode and PC camera mode cannot be used normally. This is not a product failure.
7. Do not touch the lens or other optical parts with your hands, otherwise it will cause image blur and affect imaging quality.
8. Do not disassemble the product or component structure by yourself to avoid abnormal equipment use.

### BATTERY INSTRUCTIONS

1. Please use a safe and compliant 5V/2A home charger to charge this device.
2. This device supports charging while in use.
3. If the device is left unused for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-discharging the battery.

## FUNCTION INTRODUCTION

### RIGHT SELECTION KEY

Scroll down

### MODE SWITCH KEY

Short press to switch photo/video, file browsing mode / long press to switch different color display modes

### LED LIGHT SWITCH

Left for off, right for on

### RESET BUTTON

When the product crashes, press to force shutdown

### TF CARD SLOT

Support expansion 0-64G TF card

### CHARGING INDICATOR

When charging the red light is on, when it's done the green light is on

### CHARGING INTERFACE

Input 5V=2A

### LIGHT POWER SUPPLY PORT

Output 5V=1A

### LIGHT CHARGING PORT

Input 5V=1A

### BRIGHTNESS ADJUSTMENT

Left for dim, right for bright

### RIGHT SELECTION KEY

Confirm selection / screen flip

### MENU KEY

Set the relevant parameters

### OK KEY

Confirm selection / screen flip

## WORKBENCH INTRODUCTION

## PRODUCT SPECIFICATIONS

Lens pixels	12.0MP
Magnification	1600X
Focusing mode	Manual
Fixed focus distance	90mm
Video resolution	1920*1080P
Photo resolution	4032*3024P
Video format	MOV
Photo format	JPG

## OTHER

Screen size	10.1" IPS
Screen resolution	1024*600P
Operating systems	Windows XP/ Vista/7/8/10/11, MacOS X 11.0 or higher
TF capacity	Support expansion 0-64GB
Input voltage & current	5V/2A
Battery capacity	3000mAh Lithium battery
Operating time	1 hrs
Charging time	3 hrs
Operating temperature	-10°C ~ 50°C
Camera working temperature	-15°C ~ 75°C
Light source	30 high brightness LEDs (brightness adjustable)
Languages	Simplified Chinese / Traditional Chinese / Japanese / English / Korean / French / German / Spanish / Italian / Russian / Polish

## STEPS TO USE

### INSERT THE TF CARD

Please insert the TF card when power off.

Please be careful not to mix the insertion direction and the front and back of the TF card.

Do not use excessive force when inserting or pulling out the TF card, otherwise the structure of the card slot may be damaged.

## FUNCTION BUTTON

### SETTING PAGE OPERATION INSTRUCTIONS

1. Double click the "MENU" key to enter the setting page.

2. Click "←" "→" to select the required menu item.

3. Select the required item, click the "OK" key to enter, select "←" "→" the desired parameters with and then click the "OK" key to confirm.

4. When finish setting, Click the "←" "→" key to exit the setting mode.

5. The first time using the device, you need to go to "MENU" to select the required language (the factory setting of the device is: English).

### FILE BROWSING MODE OPERATION INSTRUCTIONS

1. Click "←" twice to enter the file browsing mode, you can view all the saved photos and videos. If you need to delete the current photo or video file, click "MENU" key to enter the delete page.

2. If the current file is in video format, click "OK" key for play/pause function, click left and right key for fast forward and fast rewind function.

## MODE INTRODUCTION

The product has three mode functions, which are as follows:  
 1. Memory mode: connect to PC, read, write and format TF card.  
 2. PC camera mode: Connect to a computer and use it through computer.  
 3. Normal mode: normal use mode of the device.

### MEMORY MODE

Connect the microscope to the computer via USB data cable, select memory mode, then click the OK button to view saved photos and videos on the computer.

### PC CAMERA MODE

Only supports USB3.0 or above interface, does not support USB2.0 interface

Note: When using PC mode.

1. some function keys will not be available.
2. In PC camera mode, it is normal that the device monitor will not display images.
3. When connected, if the software has no image, please repeat plugging and unplugging the Type-c data cable or replace the cable.

## WINDOWS SYSTEM

Support Windows vista/XP/7/8/10/11 or later.  
**I. Software Installation**  
 1. For WIN 7/8 systems, download and install the "AMCAP" software.  
 2. For IN 10/11 systems, search directly for the software "Windows Camera" (Note: Please disable the default laptop camera in Windows! And you must change the camera's privacy settings, which is needed to allow access).  
 3. Or download the application "Smart Camera" for your computer, download link: [www.inskam.com/download/camera.zip](http://www.inskam.com/download/camera.zip)  
 4. Connect your microscope to the USB port on the back of your computer case with the Type-C cable that comes with the microscope, press and hold the power button to turn on your computer, select to enter PC camera mode and open the computer software to use.  
 5. Different Windows systems are recommended to use different software, please choose to download the corresponding software to use according to the system.

### MACOS SYSTEM

Support MacOS: MacOS X 11.0 or higher.  
**II. the software page display:**  
 1. Software "AMCAP" page, as follows.  
 2. "Windows Camera" page, as follows.  
 3. Software "Smart Camera" page, as follows.

## FAQ HELP

1. **The product cannot be switched on.**  
 Confirm whether the power is connected.  
 Confirm that the power supply is not sufficient due to the computer being used for charging.
2. **The product is dead.**  
 Press the device reset button to force shutdown.
3. **The image is not clear.**  
 Whether the object being measured is correctly positioned directly below the lens.  
 Whether the focal length is adjusted to the optimum position?  
 Whether the light is used correctly?
4. **Product charging.**  
 Please use 5V-2A adapter to charge the device.
5. **TF card cannot be recognized.**  
 Please make sure that the TF card format is FAT32. If not, please format it. If the device does not recognize the card, please reinsert it. The device only supports high-speed storage cards of Class 10 or above, not including Class 10.
6. **When using PC camera mode, the software does not show the picture after connecting to the computer.**  
 Reconnect the product to the computer.  
 Check if you clicked into PC mode.  
 Whether the battery is low.  
 Make sure that you are using the correct data cable and that you are not using a charging cable without data transfer function.

## COPYING PICTURE AND VIDEO FILES

### CHARGING GUIDE

Connect the device to the Type-C / USB adapter for charging (5V=2A), fast charging is not supported.

During charging, the indicator shows a red light.

When the device is fully charged, the indicator shows a green light.

When the battery is low, both memory mode and PC camera mode cannot be used normally,

318-C

**PRODUCT DESCRIPTION**  
Thank you for purchasing this equipment, it is complex to use, so please read this manual in detail to understand how to operate the microscope before using it.

**WICHTIGE HINWEISE**  
1. Please fully charge the product before first use. Please do not charge directly through the PC, but choose a 5V-2A adapter.  
2. In PC mode, this device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

**WICHTIGE HINWEISE**  
1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Power button: Press and hold for 3 seconds to switch on/off the device.  
2. Menu key: Confirm selection / screen flip.  
3. Right selection key: Confirm selection / screen flip.  
4. Left selection key: Confirm selection / screen flip.

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1. Right selection key: Confirm selection / screen flip.  
2. Left selection key: Confirm selection / screen flip.  
3. Power button: Press and hold for 3 seconds to switch on/off the device.  
4. Menu key: Confirm selection / screen flip.

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1. After opening the package, put the 3 base plates together as shown in the picture below.  
2. Combine the bracket bars and the base plates together with screws as shown in the picture below.

**FUNKTIONSTASTE**

1. After switching on, the product is in photo mode by default. Click the "OK" key to take photos, click the "Z" key once to enter the camera mode, click "ESC" to return to the main menu.  
2. Short press the "Z" key to set as: "Photographic mode" - "Video mode" - "File browsing mode", with three modes circling.  
3. In "Photography mode", click the "MENU" key to enter the menu page. In "Video mode", click the "MENU" key to enter the menu page. In "File browsing mode", click the "MENU" key to enter the menu page.

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**Notiz** If the following situation occurs for a long time during the operation, please refer to the following instructions:

a. Firstly, please disassemble the Type-C cable and reinsert it, or plug the data cable into the USB port at the back of the computer case to ensure that the device is powered properly.  
b. Secondly, if the software still cannot be installed, please check the system configuration to see if the software installation can be performed.

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**FCR-Erklärung**

Any changes or modifications not expressly approved by the party responsible for compliance may void the user's authority to operate the device.  
This device complies with Part 15 of the FCC Rules. Changes or modifications not expressly approved by the party responsible for compliance may void the user's authority to operate the device.  
This device complies with Part 15 of the FCC Rules. Changes or modifications not expressly approved by the party responsible for compliance may void the user's authority to operate the device.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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
1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

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1. Please use a clean and conforming 5V/2A helmeted adapter to charge the device.  
2. This device only supports USB3.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device cannot be used normally. This is not a product failure. If the device cannot be used normally, please contact the manufacturer for repair or replacement.  
6. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.  
7. Do not touch the lens or the camera lens with your hands. Otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or component structure yourself. Otherwise it will cause image blur and direct imaging quality.

318-C

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**PRODUCT DESCRIPTION**  
Thank you for purchasing this equipment, it is complex to use, so please read this manual in detail to understand how to operate the microscope before using it.

**REMARQUES IMPORTANTES**  
1. Please fully charge the product before first use. Pila fullness ne le fait pas charge directly through the PC, but chose a 5V-2A adapter.  
2. Dans PC mode, this device only supports USB2.0 or supérieur interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. Get appareil ne peut pas être utilisé avec précision.  
6. When using the microscope, the specific magnification effect is subject to the actual picture taken.  
7. Lorsque la batterie est faible, le mode de stockage et le mode PC camera mode cannot be used normally. This is not a product failure.  
8. Do not disassemble the product or components of an element by yourself to avoid damage to the product.

**BATTERIE INSTRUCTION**  
1. Please use a domestic 5V/2A charger to charge the device.  
2. Cet appareil prend en charge le chargement pendant son utilisation.  
3. If the device is left unused for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-charging of the battery.

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2. Combine the bracket bars and the base plates together with screws as shown in the picture below.

3. Attach the main unit to the bracket as shown in the picture below.

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After switching on, the product is in photo mode by default, click the "OK" key to take photos, click the "L" key once to enter the camera mode, click "D" for camera menu, click "M" for video mode navigation.

Short press the "L" key to set as: "Photographic mode" - "Video mode" - "File browsing mode", with three modes intuitif.

Donner l'élément de menu requis.

3. Select the required item, click the "OK" key to enter selected item, click "M" key to return to the main menu.

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

**NOTE IMPORTANT**  
1. Please fully charge the product before first use. Pila facilità no charge directly through the PC, but choose a 5V-2A adapter.  
2. The PC cable, this device only supports USB3.0 or superiore interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This device is not a precision microscope, the specific magnification effect is subject to the actual picture taken.  
6. When the battery is in charging mode and PC camera mode cannot be used normally. This is not a product failure.  
7. Do not touch the objective or other parts of the microscope, otherwise it will cause image blur and direct imaging quality.  
8. Do not disassemble the product or contact structural part of yourself to avoid damage to the microscope.

**PASTELLA ISTRUZIONI**  
1. Please use the domestic safe and conform to a 5V/2A to charge this device.  
2. Questo dispositivo supporta la ricarica durante l'uso.  
3. If the device is left unused for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-discharge of the battery.

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











318-C

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**PRODUCT DESCRIPTION**  
Thank you for purchasing this equipment. It is complex to use, so please read this manual in detail to understand how to operate the microscope before using it.

**중요 사항**  
1. Please fully charge the product before first use. Please charge directly through the PC, but choose a 5V-2A adapter.  
2. PC charging: this device only supports USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. When using the microscope, the specific magnification effect is subject to the actual picture taken.  
6. The storage mode and PC camera mode cannot be used normally.  
7. Please do not disassemble the product or connect it to other equipment, otherwise it will cause image blur and direct image quality.  
8. Do not disassemble the product or connect it to other equipment, otherwise it will cause image blur and direct image quality.

**안전 지침**  
1. Please use the product in a safe and stable environment.  
2. If the device is left unused for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-discharge of the battery.

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**Power button**  
Press and hold for 3 seconds to switch on/off the device.

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**VIKTIGA ANMÄRKNINGAR**  
1. Please fully charge the product before first use. Plug the power cable into the power socket and charge directly through the PC, but choose a 5V-2A adapter.  
2. PC include this device only supports USB2.0 or higher interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. This equipment is not a microscope, the specific magnification effect is subject to the actual picture taken.  
6. The storage mode and PC camera mode cannot be used normally. This is not a product failure.  
7. Do not touch the objective lens and the camera lens directly, otherwise it will cause image blur and direct image quality.  
8. Do not disassemble the product or component structure of yourself, ignore the safety of the equipment.

**SMETNINGAR**  
1. Please open the lid and connect the 5V/2A power adapter to the power port of the device.  
2. This device should be charged for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-discharge of the battery.

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


**Power button**  
Press and hold for 3 seconds to switch on/off the device

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**Meny nyckelvant**  
Set the parameters

**OK-knapp / Skärmen**  
Confirm selection / screen flip

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**UWAGI**  
1. Please fully charge the product before first use. Please do not charge directly through the PC, but choose a 5V-2A adapter.  
2. In computer mode, please use only supported USB 3.0 or newer interface, not USB2.0 interface.  
3. This product is not a biological microscope and is not suitable for cell biology research.  
4. This equipment is equipped with a coaxial lens. Unlike traditional digital microscopes, it maintains clarity throughout the entire magnification range without the need for refocusing. Please read the instructions carefully before use.  
5. The magnification effect is subject to the actual picture taken.  
6. When the battery level is low, the device will automatically stop working. Please recharge the battery in time.  
7. Do not disassemble the product or change the internal structure without the permission of the manufacturer. Otherwise it will cause image blur and direct damage to the device.  
8. Do not use excessive force when inserting or pulling out the TF card, otherwise it will cause damage to the card slot.

**INSTRUKCJA**  
1. Please read this manual in detail before using the product.  
2. To ensure the device is left unused for a long period of time, please ensure that it is charged at least once every 3 months to prevent irreparable damage caused by over-discharge of the battery.

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5. After fixing the workstation, press and hold the power button for 3 seconds to turn on the device and start working.

6. In Photography mode, click the "OK" key to enter the menu, click the "MENU" key to enter the "File Browsing Mode", and click the "OK" key to enter the "File Browsing Mode".

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• After switching on, the product is in photo mode by default. Click the "OK" key to take photos, click the "MENU" key to enter the camera mode, click the "OK" key to enter the camera mode, click the "MENU" key to enter the camera mode, click the "OK" key to enter the camera mode.

• Short press the "OK" key to set as: "Photographic mode" - "Video mode" - "File browsing mode", with three modes available.

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