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# Professional acoustic oral care Product user manual

## Product Description:

Sonic toothbrush is an electric toothbrush with 10,000 vibrations. The powerful motor can drive the brush to vibrate 31,000 to 40,000 times per minute. High-frequency vibration can easily remove germs and impurities around the teeth, and its cleaning range can cover all parts of the teeth. The acoustic energy produced by the motor can be transmitted to the teeth and gum surface through the brush hair. On the one hand, it can remove plaque, loosen the bond between the dental calculus and the teeth, destroy the bacteria in the gum pocket and tooth surface hidden parasite reproduction, but also can inhibit flesh inflammation and gum bleeding, prevent gingival atrophy.

## Function description:

1. Every time you finish brushing your teeth, remove the head of the toothbrush and clean it.
2. Keep the toothbrush charging interface clean and water-free when charging.
2. This product is suitable for people over 6 years old.
3. Waterproof level of this product: IP67
4. Children under 3 years of age shall not be allowed to play with this product.
5. Do not use water or acidic liquid above 40 degrees to clean the charging seat, mainframe power cord and plug, let alone immerse the plug in water.
6. Charging mode: DC direct charging, 5V = 2A
7. If your toothpaste contains peroxide, sodium bicarbonate or other bicarbonates (common ingredients in whitening toothpaste), be sure to thoroughly clean the toothbrush head with soapy water after each use. It also prevents possible plastic rupture.

## Product Description:

Sound toothbrush is an electric toothbrush with 10,000 vibrations. The powerful motor can drive the brush to vibrate 31,000 to 40,000 times per minute. High-frequency vibration can easily remove germs and impurities around the teeth, and its cleaning range can cover all parts of the teeth. The acoustic energy produced by the motor can be transmitted to the teeth and gum surface through the brush hair. On the one hand, it can remove plaque, loosen the bond between the dental calculus and the teeth, destroy the bacteria in the gum pocket and tooth surface hidden parasite reproduction, but also can inhibit flesh inflammation and gum bleeding, prevent gingival atrophy.

## Functional introduction:

1. The product has six working modes (white, clean, nurse, care, pulse, polish).
2. Switch-on within 5S is the switch of working mode, greater than 5S shutdown.
3. When the machine is working, the reminder function is every 30S, and the machine shuts down automatically for 2 minutes.
4. When the charging indicator flashes, it reminds us that the power is not needed to be charged. At this time, the power need not worry. It can also work for 10 cycles. When charging, it is necessary to clean the water at the charging interface of the electric toothbrush.

## Matters needing attention:

1. before using the product, please read the user manual carefully and keep the instructions for future reference. This product is suitable for children over 6 years of age. The following basic safety precautions should be followed when using.
2. the brush head is recommended to be replaced once every 3 months.
3. if you have had oral or gum surgery in 2 months, consult your doctor before using this product.
4. consult your dentist if you are bleeding a lot after using this product or bleeding after 1 weeks.
5. this product is only designed for cleaning teeth, gums and tongue. Do not use for any other purpose. If you feel any discomfort or pain, please stop using this product and contact your doctor.
6. Please stop using deformed or curved brush head, please change the brush head every three months, if there is wear, should be replaced in advance.
7. do not use other toothbrushes recommended by non manufacturers.
8. children should be guided to brush their teeth until they develop good brushing habits.

## Safety tips:

This product is a dry wave battery powered acoustic electric toothbrush. Battery type: USB  
bristling from DuPont, USA.

Please pay attention to the following items during use.

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