Use method and precautions of blackhead suction device:

How to use:

1. Clean your skin: Before using a blackhead blotter, clean your face with a gentle cleanser to remove oil and dirt from the surface.2. Open pores with a hot compress: Apply a hot towel to the blackhead area that needs treatment, or steam your face with a steam engine to help open pores and make blackheads easier to remove.3. Use suction function: Put the suction head of the blackhead suction device close to the skin, gently move, and use the appropriate suction to suck out the blackhead. Avoid pressing too hard or excessive suction to avoid damaging the skin.

Note:1. Before using the blackhead suction device, the skin should be cleaned first, and a hot compress can be used to open the pores and make the blackheads easier to remove. 2. Avoid sensitive areas: Do not use blackhead suction devices in sensitive areas such as eyes, eyebrows and lips, where the skin is thinner and prone to injury.3. Do not overuse: do not use the blackhead device frequently, it is generally recommended to use once or twice a week, overuse may lead to skin damage or sensitivity.How to maintain: Avoid overheating, the blackhead suction device should be placed in a well-ventilated place when charging to avoid a high temperature environment.

After charging, unplug the power in time to prevent the battery from being charged for a long time.

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