

Electronic baby Scale Model:PT-606

For the purpose of using and maintaining the baby scale, please read this instruction completely before operating

SPECIFICATIONS

Equipped with a high precision " stain-gauge " sensor and high capability CPU to ensure the weighing accuracy.

Capacity:10g-20KG Tare function
Auto power off Show low voltage

TECHNIC PARAMETER

1.Capacity:10g-20KG 2.UNIT:KG/LB/ST 3.Show numerical value:10g

4.Battery:2pcs AAA Battery 5.work temperature:0-40

6.display usage:LCD Show 3 units.

KEYFUNCTION

ON/OFF/TARE key:Power on,Power off,Power tare

UNIT key:select the measuring unit for kg/lb/st

OPERATION

Take out the scale, and open up the battery cover, install the batteries to the battery case.

Put the scale on a flat table, and ensure it steady.

WEIGHTING STEP

① Under the state press the key of " ON/OFF/TARE " to power on, the scale will be zero " 0.00 " appear on the display after 2 seconds.

② Press the UNIT key, you can change " kg ", " lb ", " st " .

③ Under the weight state, put the baby on the scale, LCD will show the weight of baby

④ Then push the unit key, to fix the weight

⑤ Press the key of " ON/OFF/TARE " , LCD shows " HELLO " then show " 0 " . Clear the weight of the container, then put the baby into container, the weight will appear on the display.

⑥ Press the key of " ON/OFF/TARE " again, power off.

ERROR DISPLAY INTRODUCTION

When the scale is over loaded, the " O-LD " will appear on the display

When the screen shows LO , it means that the power isn't enough, please change the batteries as soon as possible.

INTRODUCTION OF USING AND INSTALLING THE BATTERIES

Use 2pcs AAA Battery

When change the batteries, please open the battery cover on the back, put 2PCS AAA batteries in the battery case (make the positive poles towards), then close the battery cover tight

The product use 2PCS AAA batteries, please install it before using

ADVICE FOR USE AND CARE

Put the scale on a clean and flat table.

The max. weight of the scale is 20kg, which included the weight before tare.

To ensure the life of the scale, do not put the scale under wet , hot & cold condition.

The product belongs to high precision electronic weighing apparatus, please do not press or drop the scale otherwise the scale is easy to damage.

If the scale surface is dirty, please clean the scale with soft cloth, do not use any chemical detergent and allow water to get inside the scale.

The product is for family weighing only , do not use for legal trade.

If the result appears improperly, please check about whether the batteries and the touching point contact well.

Replace for new batteries if used for a long time already, then try again. If necessary, contact the local franchisers or telephone company.