USER MANUAL

Rechargeable
TEETH WHITENING KIT
Wireless LED

START TO WHITEN

SHELF LIFE

36 Months

HOW TO USE

1. Rinse mouthpiece and your mouth before use.
2. Looking in the mirror or cell phone to compare your teeth with shade guide paper to record which stage you are in.
3. Twist teeth whitening pen, until the gel flows to the brush.
4. Apply the gel evenly on your teeth, both upper and lower. Do not brush too much gel, more gel won't give you a faster or better result. Each time use about 0.5ml-0.7ml. Avoid applying the gel directly on your gum or inside of your lips! It will create a burning feeling. This is nothing to be alarmed about as the burning feeling will go away within 24 hours. Just make sure to rinse with water.
4. Make sure full charge before first use. Put mouthpiece into your mouth and bite it, then turn on the light. Take **deep breaths before treatment to decrease saliva production. (Excessive amount of saliva will dilute the gel's concentration.)**
5. The whitening device will turn off automatically after 16 minutes. Take it out, rinse the mouthpiece and your mouth with water and compare your teeth with shade guide paper again to see how many stages you have improved.
6. Do not eat or drink anything except water for 30 minutes after you finished the whitening.

INGREDIENTS

Propylene Glycol, Glycerin, Deionized Water, Carbomer, Hydrogen Peroxide, Cellulose Gum, Poly Vinylpyrrolidone, Menthol, Sodium Hydroxide.

STORAGE

Avoid sunlight and keep in cool place.

CAUTION

1. Not suitable for ceramic or false teeth.
2. Not suitable for teeth discoloration caused by lesion or medicine.
3. Not suitable for severe tetracycline and decayed teeth.
4. Not suitable for defective enamel, exerted dentin and damaged teeth.
5. Not suitable for children under 16 years old and pregnant women.

WARM PROMPT

1. Rinsing thoroughly, water on the teeth's surface should be dried entirely, applying gel to teeth evenly.
2. Excessive saliva will dilute the gel's concentration.
3. Avoid directly applying teeth whitening gel to gum area- wipe off excess.
4. If gel is on your face or clothes, you should wash them without any delay.
5. Please store our product in a place which is out of children's reach.

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